


Monteringsanvisning

Leggforlenger + 15 cm - Kangoo Vario



Footrest Vario incl. Enlargement mounted



Open the 4 small screws to enlarge the lower part of the footrest



Maximum enlargement by using the standard aluminium profiles

If longer adjustment is requested please order longer aluminium profiles

Take care that a longer position can clash with the small 360° turning wheel



Back view of the maximum enlargement of the footrest

Bardum

bardum.no • post@bardum.no • tlf. 64 91 8 0 60

Håndverksveien 8, 1405 Langhus