

Order Form 2025

Basic model

☒	Item no.	Size	Head circumference	Quantity
☐	C 402001	1	47 - 49 cm	
☐	C 402002	2	50 - 52 cm	
☐	C 402003	3	53 - 55 cm	
☐	C 402004	4	56 - 59 cm	
☐	C 402005	5	60 - 62 cm	

When ordering please forward "Starlight dimension sheet" fully completed.

Please select your colour

☐ grey ☐ blue ☐ bordeaux ☐ green



Net surcharge for special colour: (SB _____)

Accessories

☒	Item no.	Accessories	Quantity
☐	C 4012B	Chin protection (standard)	
☐	C 4012B	Chin protection (extra small)	
☐	C 4012B	Chin protection (extra large)	
☐	C 4013B	Forehead protection	
☐	C 4014B	Neck protection	
☐	5182-01BT	Strap system Biothane (without fastener, type of fastener freely selectable)	

Additional accessories on next page →



Starlight® Base basic model



Starlight® Base with chin, forehead and neck protection

Commission (anonymous)

Shipment: EX WORKS

Individual agreements are always superior to standard conditions.

We would like to point out, that we reserve the right to make alterations if required due to technical reasons. Possible alterations made apply also to colours and fabrics. The products displayed are either the basic models or a variation thereof (e.g. with accessories). Please, take this into consideration when ordering. We do not accept liability for printing errors. As of: August 2021

The General Terms and Conditions of ATO FORM GmbH apply.

Address/Stamp or Customer number:

Name/Company

Address

.....

.....

Customer number

Date/Signature

Order Form 2025

Accessories

☒	Item no.	Accessories	Quantity
☐	5182-07	Trial fitting (incl. minor subsequent changes, net price)	
☐	Net surcharge for express production (Next working day shipment guaranteed for orders received by 12pm.)		

Fasteners

☒	Item no.	Type	Quantity
☐	Standard	Velcro fastener	
☐	5182-02	Fixlock fastener	
☐	5182-02-SV	Fixlock fastener with safety lock	



Velcro fastener



Fixlock fastener



Fixlock fastener with safety lock

Commission (anonymous)

Shipment: EX WORKS

Individual agreements are always superior to standard conditions.

We would like to point out, that we reserve the right to make alterations if required due to technical reasons. Possible alterations made apply also to colours and fabrics. The products displayed are either the basic models or a variation thereof (e.g. with accessories). Please, take this into consideration when ordering. We do not accept liability for printing errors. As of: August 2021

The General Terms and Conditions of ATO FORM GmbH apply.

Address/Stamp or Customer number:

Name/Company

Address

.....

.....

Customer number

Date/Signature

head protection Starlight® series

Name/company

Anonymised Commission

For special dimensions and equipment, please use the:

Requirement sheet for a custom-made head protection series Starlight®

The head dimensions are determined by using a flexible measuring tape (narrow tailor tape) and are measured according to the illustrations below.




Further help on correct measuring for your head protection order can be found on our homepage at [www.https://ato-form.com/en/body-protection/head-protection-series-starlight/](https://ato-form.com/en/body-protection/head-protection-series-starlight/) or via this QR code.



General measurements

Fig. 1



*Occipital protuberance

A _____ cm

(Head circumference at eyebrow height, above the occipital protuberance)

B _____ cm

(Midway between the eyebrows to upper edge of the occipital protuberance)
Note! The B dimension is often 1 to 4 cm longer than the C dimension.

C _____ cm

(Over the highest point of the head, beginning and ending at the line of head-circumference measurement A)

D _____ cm

(Under the chin, beginning and ending at the line of head-circumference measurement A)

Important notes on filling out the form


The fitting form is only valid in conjunction with a fully completed order form.

We recommend requesting a trial fitting when ordering a head protection.

The trial fitting gives you the opportunity to make minor corrections to the measurements.

Position of the ears on the head

Fig. 2



Right side* ☐ normal

☐ low

☐ high

Left side* ☐ normal

☐ low

☐ high

*from the patient's point of view

Attention


Ideally, the tips of the ears are located on the line of measurement A.

If the positioning of the ears differs from this, make certain that measurement C is taken across the top of the head, beginning and ending at the line of measurement A.

In the case of asymmetrical ear positions, please give a clear and detailed description in the box labelled "comments" below.

Measurements position of the ears

Fig. 3



Line of measurement A to lower edge of ears

E _____ cm

(Eyebrow height (line of measurement A to lower edge of earlobes))

F _____ cm

(Distance between outer corner of the eye and measurement line E)

Ear measurement

Fig. 4

G _____ cm
(Height of the ears)


H _____ cm
(Width of the ears)



Dimension sheet for made to measure
head protection Starlight® series

To be completed if a cheek protector is required!

Fig. 5



Size and position of the
cheek protector


L _____ cm
(Height of the free area from forehead
to cheek protection)

M _____ cm
(Height of the cheek protection)

N _____ cm
(Width of the cheek protection)

Skull perspectives

Fig. 6



In the case of anatomical abnormalities,
describe the shape of the head as
precisely as possible.

If necessary, please print and alter the
forehead and/or skull shape in the
adjacent illustrations or, if applicable,
note any other deviations.

Photos are very helpful here.

Adress/Stamp or Customer number:

Name/Company

Adress

Phone

Customer number

Date/Signature